



Handheld Introduction Device Microcrystalline Pen Stainless Steel Material

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: WZ-3
- Minimum Order Quantity: 1
- Price: USD/60-100
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



Product Specification

- Name: Microcrystalline Pen
- Package Contents: Microneedle Therapy Machine, Power Adapter, User Manual
- Treatment Area: Face, Body
- Speed Levels: 5
- Function: Skin Rejuvenation, Scar Removal, Wrinkle Removal
- Power Supply: AC 100-240V, 50/60Hz
- Usage: Home, Salon
- Needle Length: 0.25mm-2.5mm
- Power Consumption: 10W
- Highlight: Handheld Microcrystalline Pen, Introduction Device Microcrystalline Pen



More Images



Product Parameters



Specification

| | | | |
|---------------------------|--|-----------------|--|
| Brand | Hydra pen | Model | H5 |
| EMS intensity level | 4 levels | EMS frequency | 0.25 Hz |
| Loadable liquid capacity | 3 ml | Speed Level | 6 Speeds |
| Adjustable liquid output | 3 levels | Handle Material | ABS |
| Needle Adjustment Range | 0-1.0 mm | Adapter | US/EU/UK/AU/KR |
| Battery capacity | 700mAh | G.W. | 19.4*14.6*6.9 cm/1 pc/ 690 g |
| Available Catridges | H12, H36, H24, Nano-HS, Nano-HR | 1 carton | 41*32*38 cm/20 pcs/? kg |
| Customization (OEM & ODM) | logo, private label, package, new shape, change color. | Packing List | ① Box ② Inductive base ③ Hydra pen H5 ④ Cable ⑤ Adapter ⑥ Manual ⑦ Retaining ring* 2 pcs ⑧ 12 pins * 2 pcs |

4-pole release micro current

Using 4-pole micro current to release molecules that penetrate the bottom layer.



6-speed deep penetration

1-6 speed intelligent speed regulation solves various skin problems.



LED display screen intelligent pause

LED digital display can be used for charging and plugging, making operation more convenient.



Adjustable needle length 0~1.0 mm

Deep adjustment, The functions of regulation, repair, and hydration are fully equipped.





EMS | infusion carrier essence device





Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou